

CD58 Mouse anti-Human Monoclonal (RPE) (TS2/9) Antibody - LS-C134513 - LSBio

<b>CatalogID:</b>	LS-C134513
<b>Target:</b>	CD58 molecule
<b>Synonyms:</b>	CD58 Antibody, CD58 antigen Antibody, AG3 Antibody, LFA-3 Antibody, LFA3 Antibody, Surface glycoprotein LFA-3 Antibody, Surface glycoprotein lfa3 Antibody, CD58 molecule Antibody
<b>Family / Subfamily:</b>	Immunoglobulin / not assigned-Immunoglobulin
<b>Host</b>	CD58 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	TS2/9
<b>Conjugations:</b>	R. Phycoerythrin (RPE)
<b>Immunogen Species:</b>	CD58 antibody was raised against Human
<b>Immunogen:</b>	CD58 antibody was raised against human cytotoxic T cells.
<b>Specificity:</b>	Both transmembrane and GPI anchored CD58.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Gel filtration
<b>Presentation:</b>	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
<b>Recommended Storage:</b>	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
<b>Usage Summary:</b>	Function: Inhibits HLA-DR mediated T cell cytotoxicity.
<b>Uses:</b>	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
<b>Size:</b>	120 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences