


CD63 Mouse anti-Human Monoclonal (RPE) (AHN16.1/46-4-5) Antibody - LS-C134527 - LSBio

CatalogID:	LS-C134527
Target:	CD63 molecule
Synonyms:	CD63 Antibody, CD63 molecule Antibody, LAMP-3 Antibody, ME491 Antibody, Granulophysin Antibody, MLA1 Antibody, Tetraspanin-30 Antibody, Tspan-30 Antibody, CD63 antigen Antibody, OMA81H Antibody, TSPAN30 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD63 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	AHN16.1/46-4-5
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD63 antibody was raised against Human
Immunogen:	CD63 antibody was raised against human eosinophils.
Specificity:	CD63 molecule of about 50 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	ICC, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst

Flow Cytometry Image:



Flow cytometry of CD63 antibody

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences