

IFN Gamma / Interferon Gamma Mouse anti-Human Monoclonal (RPE) (ANC2E11) Antibody - LS-C134584 - LSBio	
CatalogID:	LS-C134584
Target:	interferon, gamma (IFNG)
Synonyms:	IFNG Antibody, IFG Antibody, Immune interferon Antibody, IFI Antibody, Interferon gamma Antibody, IFN-gamma Antibody, Interferon, gamma Antibody
Family / Subfamily:	Interferon / not assigned-Interferon
Host	IFNG antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	ANC2E11
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	IFN Gamma / Interferon Gamma antibody was raised against Human
Immunogen:	IFN Gamma / Interferon Gamma antibody was raised against recombinant human IFN-gamma.
Specificity:	Human interferon gamma.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Flow Cytometry: Intracellular staining.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/25/2014 © 2014 LifeSpan BioSciences	