

SELL / L-Selectin / CD62L Mouse anti-Human Monoclonal (FITC) (LAM1-116) Antibody - LS-C134600 - LSBio	
CatalogID:	LS-C134600
Target:	selectin L (SELL)
Synonyms:	SELL Antibody, CD62L Antibody, CD62L antigen Antibody, LAM-1 Antibody, LAM1 Antibody, LECAM1 Antibody, LEU8 Antibody, Leukocyte adhesion molecule 1 Antibody, Lyam-1 Antibody, LSEL Antibody, gp90-MEL Antibody, HLHRc Antibody, PLNHR Antibody, L-selectin Antibody, Leu-8 Antibody, Selectin L Antibody, LNHR Antibody, LYAM1 Antibody, Lymph node homing receptor Antibody, Lymphocyte adhesion molecule 1 Antibody, TQ1 Antibody, PIn homing receptor Antibody
Host	SELL antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	LAM1-116
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SELL / L-Selectin / CD62L antibody was raised against Human
Antigen Type:	Cells
Immunogen:	SELL / L-Selectin / CD62L antibody was raised against mouse pre-B cells expressing human L-selectin.
Specificity:	Lectin binding domain of human, mouse, rat and guinea pig CD62L.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Blocks CD62L function and induces expression of beta-1 and beta-2 integrins.
Uses:	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	