


TNFRSF18 / GITR Mouse anti-Human Monoclonal (FITC) (ANC7D6-1B7) Antibody - LS-C134624 - LSBio	
CatalogID:	LS-C134624
Target:	tumor necrosis factor receptor superfamily, member 18 (TNFRSF18)
Synonyms:	TNFRSF18 Antibody, AITR Antibody, CD357 Antibody, CD357 antigen Antibody, GITR Antibody, GITR-D Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	TNFRSF18 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgM,k
Clone Name:	ANC7D6-1B7
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TNFRSF18 / GITR antibody was raised against Human
Immunogen:	TNFRSF18 / GITR antibody was raised against recombinant human GITR, PHA stimulated human PBL.
Specificity:	Human GITR in EIA and on cell surface in FACS. Preincubation of the receptor with an excess of this mAb will block binding of recombinant GITR Ligand.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	120 tst

Flow Cytometry Image:



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences