


TNFRSF14 / CD270 / HVEM Mouse anti-Human Monoclonal (FITC) (ANC3B7) Antibody - LS-C134628 - LSBio

CatalogID:	LS-C134628
Target:	tumor necrosis factor receptor superfamily, member 14 (TNFRSF14)
Synonyms:	TNFRSF14 Antibody, CD270 Antibody, CD270 antigen Antibody, CD40-like protein Antibody, Herpesvirus entry mediator Antibody, Herpesvirus entry mediator A Antibody, HVEM Antibody, LIGHTR Antibody, Herpes virus entry mediator A Antibody, TR2 Antibody, ATAR Antibody, HVEA Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	TNFRSF14 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	ANC3B7
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TNFRSF14 / CD270 / HVEM antibody was raised against Human
Immunogen:	TNFRSF14 / CD270 / HVEM antibody was raised against recombinant human CD270(HVEM).
Specificity:	Recombinant CD270 in EIA, and binds to CD270 present on human U-937 cell surface in Flow Cytometry.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	120 tst

Flow Cytometry Image:



Flow cytometry of CD270 / TNFRSF14 / HVEM antibody

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences