

CD40 Mouse anti-Human Monoclonal (RPE) (EA-5) Antibody - LS-C134642 - LSBio

CatalogID:	LS-C134642
Target:	CD40 molecule, TNF receptor superfamily member 5
Synonyms:	CD40 Antibody, B-cell surface antigen CD40 Antibody, B cell-associated molecule Antibody, B cell surface antigen CD40 Antibody, CD40 antigen Antibody, CD40 type II isoform Antibody, p50 Antibody, Bp50 Antibody, CD40L receptor Antibody, CDW40 Antibody, TNFRSF5 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	CD40 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	EA-5
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD40 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD40 antibody was raised against murine L cells transfectants expressing human CD40.
Specificity:	CD40 molecule of about 44 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Binding of cross linked EA-5 will partially activate B cells.
Uses:	Flow Cytometry, ELISA, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences