

CD84 Mouse anti-Human Monoclonal (FITC) (152.1D5) Antibody - LS-C134661 - LSBio

CatalogID:	LS-C134661
Target:	CD84 molecule
Synonyms:	CD84 Antibody, CD84 antigen Antibody, Cell surface antigen MAX.3 Antibody, CD84 molecule Antibody, HCD84 Antibody, Hly9-beta Antibody, LY9B Antibody, Leukocyte antigen CD84 Antibody, SLAMF5 Antibody, MCD84 Antibody, SLAM family member 5 Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD84 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	152.1D5
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD84 antibody was raised against Human
Immunogen:	CD84 antibody was raised against human MLR stimulated lymphocytes.
Specificity:	Molecule of about 74 kd.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences