

CD86 Mouse anti-Human Monoclonal (FITC) (BU63) Antibody - LS-C134666 - LSBio

CatalogID:	LS-C134666
Target:	CD86 molecule
Synonyms:	CD86 Antibody, Activation B7-2 antigen Antibody, B7-2 Antibody, B7.2 Antibody, B70 Antibody, CD86 molecule Antibody, CTLA-4 counter-receptor B7.2 Antibody, CD86 antigen Antibody, FUN-1 Antibody, LAB72 Antibody, BU63 Antibody, CD28LG2 Antibody, Ctla-4 counter-receptors b7.2 Antibody
Host	CD86 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	BU63
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD86 antibody was raised against Human
Immunogen:	CD86 antibody was raised against aRH-77 (B-lymphoblastoid cell line).
Specificity:	CD86 molecule of 80 kd.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Blocks MLR and blocks binding of soluble CD152-mulg fusion protein to CD86.
Uses:	Western blot, Flow Cytometry, ELISA, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences