

KLRD1 / CD94 Mouse anti-Human Monoclonal (FITC) (HP-3D9) Antibody - LS-C134677 - LSBio	
CatalogID:	LS-C134677
Target:	killer cell lectin-like receptor subfamily D, member 1 (KLRD1)
Synonyms:	KLRD1 Antibody, CD94 Antibody, KP43 Antibody, NK cell receptor Antibody, CD94 antigen Antibody
Host	KLRD1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	HP-3D9
Conjugations:	Fluorescein (FITC)
Immunogen Species:	KLRD1 / CD94 antibody was raised against Human
Immunogen:	KLRD1 / CD94 antibody was raised against cultured human NK cells.
Specificity:	CD94 molecule of 70 kd.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Inhibits IL-2 dependent proliferation of NK cells.
Uses:	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014 © 2014 LifeSpan BioSciences	