

CD98 Mouse anti-Human Monoclonal (FITC) (UM7F8) Antibody - LS-C134687 - LSBio

CatalogID:	LS-C134687
Target:	CD98
Host	CD98 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	UM7F8
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD98 antibody was raised against Human
Immunogen:	CD98 antibody was raised against molt 13 T cell line.
Specificity:	Heavy chain of a 125 kd heterodimeric membrane glycoprotein CD98.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Comitogenic with soluble anti-CD2 and immobilized anti-CD3 mAbs.
Uses:	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences