

C-Kit / CD117 Mouse anti-Human Monoclonal (RPE) (57A5) Antibody - LS-C134711 - LSBio	
CatalogID:	LS-C134711
Target:	v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog
Synonyms:	KIT Antibody, C-Kit Antibody, CD117 Antibody, PBT Antibody, Piebald trait Antibody, Piebald trait protein Antibody, Proto-oncogene c-Kit Antibody, SCFR Antibody, Soluble KIT variant 1 Antibody, Tyrosine-protein kinase Kit Antibody, p145 c-kit Antibody, CD117 antigen Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	KIT antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	57A5
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	C-Kit / CD117 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	C-Kit / CD117 antibody was raised against human MO7e tumor cells.
Specificity:	Epitope 2 of the CD117 molecule of about 150 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Does not block binding of SCF.
Uses:	Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	