

BCL2 / Bcl-2 Mouse anti-Human Monoclonal (RPE) (Bcl-2/100) Antibody - LS-C134745 - LSBio	
CatalogID:	LS-C134745
Target:	B-cell CLL/lymphoma 2 (BCL2)
Synonyms:	BCL2 Antibody, Apoptosis regulator Bcl-2 Antibody, B-cell CLL/lymphoma 2 Antibody, B-cell lymphoma protein 2 Antibody, Bcl-2 Antibody, PPP1R50 Antibody
Family / Subfamily:	Apoptosis
Host	BCL2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	Bcl-2/100
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	BCL2 / Bcl-2 antibody was raised against Human
Immunogen:	BCL2 / Bcl-2 antibody was raised against human bcl-2 peptide 41-54.
Specificity:	Intracellular 26 kd bcl-2 protein.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Uses:	IHC - Paraffin, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	