

TNFRSF9 / 4-1BB / CD137 Mouse anti-Human Monoclonal (FITC) (4B4-1) Antibody - LS-C134758 - LSBio	
CatalogID:	LS-C134758
Target:	tumor necrosis factor receptor superfamily, member 9 (TNFRSF9)
Synonyms:	TNFRSF9 Antibody, 4-1BB ligand receptor Antibody, 4-1BB Antibody, CD137 antigen Antibody, CD137 Antibody, Homolog of mouse 4-1BB Antibody, ILA Antibody, Receptor protein 4-1BB Antibody, T cell antigen ILA Antibody, T-cell antigen 4-1BB homolog Antibody, T-cell antigen ILA Antibody, CDw137 Antibody
Family / Subfamily:	TNF Receptor / not assigned-TNF Receptor
Host	TNFRSF9 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	4B4-1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TNFRSF9 / 4-1BB / CD137 antibody was raised against Human
Immunogen:	TNFRSF9 / 4-1BB / CD137 antibody was raised against ectodomain of recombinant human 4-1BB fusion protein.
Specificity:	CD137 molecule of about 35 kd from activated T cells.
Reactivity:	Human
Purification:	Protein A purified
Presentation:	50 mM sodium phosphate, pH 7.5, 100 mM potassium chloride, 150 mM sodium chloride, 5% Glycerol, 0.2% BSA, 0.04% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Function: Blocks binding of CD137 hulg to Raji cells. Immobilized mAb induces TNF a production by monocytes.
Uses:	Flow Cytometry, ELISA, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	120 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	