

**CD45RA Mouse anti-Human Monoclonal (RPE) (158.4D3) Antibody - LS-C134797 - LSBio**

<b>CatalogID:</b>	LS-C134797
<b>Target:</b>	CD45RA
<b>Host</b>	CD45RA antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	158.4D3
<b>Conjugations:</b>	R. Phycoerythrin (RPE)
<b>Immunogen Species:</b>	CD45RA antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD45RA antibody was raised against stimulated human leukocytes.
<b>Specificity:</b>	All CD45 molecules containing the A region exon.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Gel filtration
<b>Presentation:</b>	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 15% Glycerol, 0.2% BSA, 0.04% sodium azide
<b>Recommended Storage:</b>	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
<b>Uses:</b>	IHC - Paraffin, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	120 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences