

CD45.1 Mouse anti-Mouse Monoclonal (PE, Cy7) (A20) Antibody - LS-C13596 - LSBio

CatalogID:	LS-C13596
Target:	CD45.1
Host	CD45.1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	A20
Conjugations:	Phycoerythrin (PE), Cy7
Immunogen Species:	CD45.1 antibody was raised against Mouse
Specificity:	Recognizes mouse CD45.1 (Leukocyte Common Antigen). Recognizes all mouse leukocytes expressing the CD 45.1 allotype.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry and fluorescence microscopy. Markers for origin of donor and host cells. Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences