

CD45.2 Mouse anti-Mouse Monoclonal (PE, Cy5.5) (104-2) Antibody - LS-C13609 - LSBio

CatalogID:	LS-C13609
Target:	CD45.2
Host	CD45.2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	104-2
Conjugations:	Phycoerythrin (PE), Cy5.5
Immunogen Species:	CD45.2 antibody was raised against Mouse
Specificity:	Recognizes mouse CD45.2 (Leukocyte Common Antigen). Recognizes CD45 (Leukocyte Common Antigen) on all leukocytes of mouse strains expressing the CD45.2 allotype (e. g. A, AKR, BALB/c, CBA/Ca, CBA/J, C3H/He, C57BL, C57BR, C57L, C58, DBA/1, DBA/2, NZB, SWR, 129).
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry and fluorescence microscopy. Markers for origin of donor and host cells. Immunoprecipitation. Flow Cytometry: 0.3 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences