

CD45.2 Mouse anti-Mouse Monoclonal (FITC) (104-2) Antibody - LS-C13612 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C13612 |
| Target: | CD45.2 |
| Host | CD45.2 antibody was produced in Mouse |
| Clonality: | Monoclonal |
| Isotype: | IgG2a,k |
| Clone Name: | 104-2 |
| Conjugations: | Fluorescein (FITC) |
| Immunogen Species: | CD45.2 antibody was raised against Mouse |
| Specificity: | Recognizes mouse CD45.2 (Leukocyte Common Antigen). Recognizes CD45 (Leukocyte Common Antigen) on all leukocytes of mouse strains expressing the CD45.2 allotype (e. g. A, AKR, BALB/c, CBA/Ca, CBA/J, C3H/He, C57BL, C57BR, C57L, C58, DBA/1, DBA/2, NZB, SWR, 129). |
| Reactivity: | Mouse |
| Purification: | Purified |
| Presentation: | PBS, sodium azide. |
| Recommended Storage: | Long term: Add glycerol (40-50%) -20°C; Short term: +4°C |
| Usage Summary: | Flow cytometry and fluorescence microscopy. Markers for origin of donor and host cells. Immunoprecipitation. Flow Cytometry: 1 ug/1E6 cells. |
| Uses: | Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher) |
| Size: | 500 µg |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences