

MYC / c-Myc Rabbit anti-Human Polyclonal Antibody - LS-C136226 - LSBio

CatalogID:	LS-C136226
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	MYC / c-Myc antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MYC / c-Myc antibody was raised against synthetic peptide {EQKLISEEDL (410-419)} from human c-Myc (c Myc tag) conjugated to an immunogenic carrier protein was used as the antigen.
Specificity:	Specific for c-Myc.
Reactivity:	Human
Purification:	Purified
Reconstitution:	Sterile water
Presentation:	Lyophilized. Centrifuge to remove any insoluble material
Recommended Storage:	Maintain the lyophilized/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Uses:	Western blot (10 - 50 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences