


MYC / c-Myc Rabbit anti-Human Polyclonal Antibody - LS-C137025 - LSBio

CatalogID:	LS-C137025
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	MYC / c-Myc antibody was raised against Human
Antigen Type:	Synthetic peptide - KLH conjugated
Immunogen:	MYC / c-Myc antibody was raised against synthetic peptide conjugated to KLH
Reactivity:	Human, Mouse, Rat
Purification:	Affinity purified
Presentation:	PBS, pH 7.4, 0.02% sodium azide
Recommended Storage:	Store at -20°C for up to one year.
Uses:	IHC (1:50 - 1:200), Western blot (1:250 - 1:1000), ELISA (1:20000 - 1:80000) (Optimal dilution to be determined by the researcher)
Size:	100 µl

Western Blot Image:



Western Blot analysis of HeLa cell lysates with anti-c-Myc antibody

Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	