

ITGA5/Integrin Alpha 5/CD49e Rat anti-Mouse Monoclonal (PE) (MFR5) Antibody - LS-C13707 - LSBio	
CatalogID:	LS-C13707
Target:	integrin, alpha 5 (fibronectin receptor, alpha polypeptide) (ITGA5)
Synonyms:	ITGA5 Antibody, FNRA Antibody, Integrin alpha-5 Antibody, Integrin Alpha 5 Antibody, Integrin alpha-F Antibody, VLA-5 Antibody, CD49e Antibody, CD49e antigen Antibody, VLA5A Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGA5 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	MFR5
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ITGA5/Integrin Alpha 5/CD49e antibody was raised against Mouse
Immunogen:	ITGA5/Integrin Alpha 5/CD49e antibody was raised against c57BL/6 x A/J F1 mouse mast cell line MC/9.
Specificity:	Recognizes the alpha5 subunit of the integrin a5b1 fibronectin receptor. Reacts with the alpha5 chain (VLA-5) of the CD49e/CD29 heterodimeric fibronectin receptor.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. In vitro blocking of adhesion. Immunohistochemistry (frozen sections). T-cell co-stimulation. Flow Cytometry: 0.3 ug/1E6 cells.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/25/2014 © 2014 LifeSpan BioSciences	