

ITGA6/Integrin Alpha 6/CD49f Mouse anti-Human Monoclonal (PE) Antibody - LS-C13712 - LSBio	
<b>CatalogID:</b>	LS-C13712
<b>Target:</b>	integrin, alpha 6 (ITGA6)
<b>Synonyms:</b>	ITGA6 Antibody, CD49f Antibody, CD49f antigen Antibody, Integrin alpha6B Antibody, Integrin, alpha 6 Antibody, Integrin alpha-6 precursor Antibody, VLA-6 Antibody, Integrin alpha-6 Antibody, ITGA6B Antibody
<b>Family / Subfamily:</b>	Integrin / not assigned-Integrin
<b>Host</b>	ITGA6 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	ITGA6/Integrin Alpha 6/CD49f antibody was raised against Human
<b>Specificity:</b>	Recognizes human CD49f.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Recommended Storage:</b>	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
<b>Usage Summary:</b>	Flow cytometry. Immunoprecipitation. Immunohistochemistry (frozen sections). Flow Cytometry: 10 ul/1E6 cells.
<b>Uses:</b>	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> <b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b> Created on 9/25/2014 © 2014 LifeSpan BioSciences	