


AP1G1 / Gamma Adaptin Rabbit anti-Human Polyclonal Antibody - LS-C137137 - LSBio

<b>CatalogID:</b>	LS-C137137
<b>Target:</b>	adaptor-related protein complex 1, gamma 1 subunit (AP1G1)
<b>Synonyms:</b>	AP1G1 Antibody, ADTG Antibody, AP-1 complex subunit gamma-1 Antibody, CLAPG1 Antibody, Gamma adaptin Antibody, Gamma1-adaptin Antibody
<b>Host</b>	AP1G1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Immunogen Species:</b>	AP1G1 / Gamma Adaptin antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	AP1G1 / Gamma Adaptin antibody was raised against a peptide from human AP1G1 protein
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.4, 0.02% sodium azide
<b>Recommended Storage:</b>	Store at -20°C for up to one year.
<b>Uses:</b>	IHC (1:100 - 1:500), Western blot (1:200 - 1:1000), ELISA (1:1000 - 1:10000) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl

Western Blot Image:



<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	