


POU5F1 / OCT4 Mouse anti-Human Monoclonal (aa102-265) Antibody - LS-C137531 - LSBio

CatalogID:	LS-C137531
Target:	POU class 5 homeobox 1 (POU5F1)
Synonyms:	POU5F1 Antibody, Oct-4 Antibody, OCT3 Antibody, OCT4 Antibody, OTF4 Antibody, POU class 5 homeobox 1 Antibody, Octamer-binding protein 3 Antibody, Octamer-binding protein 4 Antibody, OTF-3 Antibody, Oct-3 Antibody, OTF3 Antibody
Host	POU5F1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Immunogen Species:	POU5F1 / OCT4 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	POU5F1 / OCT4 antibody was raised against recombinant protein corresponding to fragment (amino acids 102-265) POU5F1
Specificity:	Anti-POU5F1 recognizes recombinant protein POU5F1.
Epitope:	aa102-265
Reactivity:	Human
Purification:	Affinity purified
Presentation:	Citrate-Tris-HCl buffer, pH 7.0, 0.02% Proclin 300
Recommended Storage:	Store at -20°C for up to one year.
Uses:	Western blot (1:1000), ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl

Western Blot Image:



Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	