


NDRG2 Mouse anti-Human Monoclonal (1D12) Antibody - LS-C139347 - LSBio

CatalogID:	LS-C139347
Target:	NDRG family member 2 (NDRG2)
Synonyms:	NDRG2 Antibody, KIAA1248 Antibody, N-myc downstream regulator 2 Antibody, NDRG family member 2 Antibody, SYLD Antibody, NDR1-related protein NDR2 Antibody, Protein Syld709613 Antibody, Syld709613 protein Antibody, Cytoplasmic protein Ndr1 Antibody, Protein NDRG2 Antibody
Host	NDRG2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	1D12
Immunogen Species:	NDRG2 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	NDRG2 antibody was raised against nDRG2 (NP_057334, 1 a.a. ~ 97 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Specificity:	Human NDRG2
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2
Recommended Storage:	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.
Usage Summary:	Western Blot: Cell lysate, Recombinant protein.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg

Western Blot Image:



NDRG2 monoclonal antibody (M06), clone 1D12. Western blot of NDRG2 expression in A-431.

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences