

CBLL1 / HAKAI Mouse anti-Human Monoclonal (3B12) Antibody - LS-C139471 - LSBio

CatalogID:	LS-C139471
Target:	Cbl proto-oncogene-like 1, E3 ubiquitin protein ligase (CBLL1)
Synonyms:	CBLL1 Antibody, C-Cbl-like protein 1 Antibody, E-cadherin binding protein E7 Antibody, HAKAI Antibody, RING finger protein 188 Antibody, RNF188 Antibody
Host	CBLL1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	3B12
Immunogen Species:	CBLL1 / HAKAI antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	CBLL1 / HAKAI antibody was raised against cBLL1 (NP_079090, 1 a.a. ~ 101 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Specificity:	Human CBLL1
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2
Recommended Storage:	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.
Usage Summary:	Western Blot: Recombinant protein.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences