

SELL / L-Selectin / CD62L Rat anti-Mouse Monoclonal (PE, Cy5) (MEL-14) Antibody - LS-C13955 - LSBio	
CatalogID:	LS-C13955
Target:	selectin L (SELL)
Synonyms:	SELL Antibody, CD62L Antibody, CD62L antigen Antibody, LAM-1 Antibody, LAM1 Antibody, LECAM1 Antibody, LEU8 Antibody, Leukocyte adhesion molecule 1 Antibody, Lyam-1 Antibody, LSEL Antibody, gp90-MEL Antibody, HLHRc Antibody, PLNHR Antibody, L-selectin Antibody, Leu-8 Antibody, Selectin L Antibody, LNHR Antibody, LYAM1 Antibody, Lymph node homing receptor Antibody, Lymphocyte adhesion molecule 1 Antibody, TQ1 Antibody, PIn homing receptor Antibody
Host	SELL antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	MEL-14
Conjugations:	Phycoerythrin (PE), Cy5
Immunogen Species:	SELL / L-Selectin / CD62L antibody was raised against Mouse
Immunogen:	SELL / L-Selectin / CD62L antibody was raised against c3H/eb cloned mouse B lymphoma 38C-13.
Specificity:	Recognizes CD62L/L-selectin, Mr 90kD.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. In vivo and in vitro blocking of adhesion. Cytotoxicity studies. Flow Cytometry: 0.1 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	