

TLR2 Mouse anti-Mouse Monoclonal (FITC) (T2.5) Antibody - LS-C139893 - LSBio

CatalogID:	LS-C139893
Target:	toll-like receptor 2 (TLR2)
Synonyms:	TLR2 Antibody, CD282 antigen Antibody, CD282 Antibody, TIL4 Antibody, Toll-like receptor 2 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	T2.5
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR2 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	TLR2 antibody was raised against mouse TLR2 peptide
Specificity:	Mouse TLR2
Reactivity:	Mouse, Human
Purification:	Purified
Presentation:	PBS, 0.02% sodium azide, 1% BSA.
Recommended Storage:	Store at 4°C, stable for one year.
Uses:	IHC - Frozen, Immunofluorescence, Immunoprecipitation, Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences