

TLR2 Mouse anti-Mouse Monoclonal (FITC) (mT2.7) Antibody - LS-C139898 - LSBio

CatalogID:	LS-C139898
Target:	toll-like receptor 2 (TLR2)
Synonyms:	TLR2 Antibody, CD282 antigen Antibody, CD282 Antibody, TIL4 Antibody, Toll-like receptor 2 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	mT2.7
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR2 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	TLR2 antibody was raised against mouse TLR2 peptide
Specificity:	Mouse TLR2
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.02% sodium azide, 1% BSA
Recommended Storage:	Store at 4°C, stable for one year.
Usage Summary:	For immunohistology and flow cytometry, dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences