

TNF Alpha Mouse anti-Human Monoclonal (FITC) (52B83) Antibody - LS-C139943 - LSBio	
CatalogID:	LS-C139943
Target:	tumor necrosis factor (TNF)
Synonyms:	TNF Antibody, APC1 protein Antibody, Cachectin Antibody, DIF Antibody, TNF Alpha Antibody, Tumor necrosis factor-alpha Antibody, TNF-alpha Antibody, TNF-a Antibody, TNFA Antibody, Tumor necrosis factor Antibody, TNF, macrophage-derived Antibody, TNF, monocyte-derived Antibody, TNFSF2 Antibody
Family / Subfamily:	TNF / not assigned-TNF
Host	TNF antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	52B83
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TNF Alpha antibody was raised against Human
Specificity:	Human TNF
Reactivity:	Human, Monkey, Mouse, Guinea pig, Pig
Purification:	Purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	Store at 4°C, stable for one year.
Usage Summary:	For immunohistology, flow cytometry and Western blotting dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50.
Uses:	IHC - Paraffin, IHC - Frozen, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	