

TLR4 Mouse anti-Human Monoclonal (FITC) (HTA125) Antibody - LS-C140002 - LSBio

CatalogID:	LS-C140002
Target:	toll-like receptor 4 (TLR4)
Synonyms:	TLR4 Antibody, CD284 Antibody, CD284 antigen Antibody, Homolog of Drosophila toll Antibody, TLR-4 Antibody, TOLL Antibody, ARMD10 Antibody, Toll-like receptor 4 Antibody, HToll Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	HTA125
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR4 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	TLR4 antibody was raised against ba/F3 cells expressing TLR4
Specificity:	Human TLR4
Reactivity:	Human, Dog
Purification:	Purified
Presentation:	PBS, 0.02% sodium azide, 1% BSA
Recommended Storage:	Store at 4°C, stable for one year.
Uses:	Immunofluorescence, Western blot, Immunoprecipitation, Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences