

TLR2 Mouse anti-Human Monoclonal (FITC) (TLR2.45) Antibody - LS-C140128 - LSBio

CatalogID:	LS-C140128
Target:	toll-like receptor 2 (TLR2)
Synonyms:	TLR2 Antibody, CD282 antigen Antibody, CD282 Antibody, TIL4 Antibody, Toll-like receptor 2 Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	TLR2.45
Conjugations:	Fluorescein (FITC)
Immunogen Species:	TLR2 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	TLR2 antibody was raised against ba/F3 cells that stably expressed human Flag-tagged TLR2
Specificity:	Human TLR2
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	Store at 4°C, stable for one year.
Usage Summary:	For flow cytometry and immunofluorescence dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50.
Uses:	Immunofluorescence (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences