

ASGR1 / ASGPR Mouse anti-Rat Monoclonal (FITC) (8D7) Antibody - LS-C140170 - LSBio

CatalogID:	LS-C140170
Target:	asialoglycoprotein receptor 1 (ASGR1)
Synonyms:	ASGR1 Antibody, ASGPR1 Antibody, ASGPR Antibody, ASGPR 1 Antibody, Asialoglycoprotein receptor 1 Antibody, Hepatic lectin H1 Antibody, HL-1 Antibody, ASGP-R 1 Antibody, CLEC4H1 Antibody
Host	ASGR1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	8D7
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ASGR1 / ASGPR antibody was raised against Rat
Immunogen:	ASGR1 / ASGPR antibody was raised against crude rat liver membrane extracts
Specificity:	Rat ASGR1
Reactivity:	Rat, Human
Purification:	Purified
Presentation:	PBS, 1% BSA, 0.02% sodium azide
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	For immunohistology, flow cytometry and Western blotting dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50.
Uses:	IHC - Frozen, Immunofluorescence, Western blot, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences