

**CR1L / CRRY Mouse anti-Rat Monoclonal (FITC) (TLD-1C11) Antibody - LS-C140183 - LSBio**

<b>CatalogID:</b>	LS-C140183
<b>Target:</b>	Complement Component (3b/4b) Receptor 1-like (CR1L)
<b>Synonyms:</b>	CR1L Antibody, Cr1l Antibody, Crry Antibody
<b>Host</b>	CR1L antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	TLD-1C11
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	CR1L / CRRY antibody was raised against Rat
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CR1L / CRRY antibody was raised against rat microglial cells
<b>Specificity:</b>	Rat
<b>Reactivity:</b>	Rat
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% BSA
<b>Recommended Storage:</b>	Store at 4°C, stable for one year.
<b>Uses:</b>	IHC - Frozen, Immunofluorescence, Western blot, Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences