

CD33 Mouse anti-Human Monoclonal (PE) (FOS) Antibody - LS-C140337 - LSBio	
CatalogID:	LS-C140337
Target:	CD33 molecule
Synonyms:	CD33 Antibody, CD33 antigen Antibody, CD33 molecule Antibody, gp67 Antibody, SIGLEC-3 Antibody, p67 Antibody, CD33 antigen (gp67) Antibody, SIGLEC3 Antibody
Family / Subfamily:	Immunoglobulin / not assigned-Immunoglobulin
Host	CD33 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	FOS
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD33 antibody was raised against Human
Immunogen:	CD33 antibody was raised against fMY9S5 cells containing the CD33 gene.
Specificity:	Human CD33
Reactivity:	Human
Purification:	Protein A/G purified
Presentation:	PBS, 0.08% sodium azide, 0.2% carrier protein, sterile-filtered
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	PBMC: Add 10 ul of antibody/10 ⁶ PBMC in 100 ul PBS. Mix gently and incubate for 15 minutes at 2 to 8°C. Wash twice with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. WHOLE BLOOD: Add 10 ul of antibody/100 ul of whole blood. Mix gently and incubate for 15 minutes at room temperature 20°C. Lyse the whole blood. Wash once with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. See instrument manufacturers instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	