

CD45RA Mouse anti-Human Monoclonal (PE) (BHKG-2) Antibody - LS-C140362 - LSBio

CatalogID:	LS-C140362
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	BHKG-2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD45RA antibody was raised against Human
Immunogen:	CD45RA antibody was raised against peripheral blood monocytes from rheumatoid arthritis patient.
Specificity:	Human PTPRC / CD45
Reactivity:	Human
Purification:	Protein A/G purified
Presentation:	PBS, 0.08% sodium azide, 0.2% carrier protein, sterile-filtered
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	Add 10 ul of antibody / 100 ul of 10 ⁶ PBMC in PBS. Mix Wash twice with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. Add 10 ul of antibody / 100 ul of Whole Blood. Mix gently and incubate for 15 minutes at room temperature 20°C. Lyse Whole Blood and analyze or fix with 0.5% v/v paraformaldehyde in PBS and analyze. See instrument manufacturers instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences