

ITGAM / CD11b Mouse anti-Human Monoclonal (FITC) (F6.2) Antibody - LS-C140417 - LSBio	
CatalogID:	LS-C140417
Target:	integrin, alpha M (complement component 3 receptor 3 subunit) (ITGAM)
Synonyms:	ITGAM Antibody, Antigen CD11b (p170) Antibody, CD11B Antibody, CD11b antigen Antibody, CR3A Antibody, Integrin alpha-M Antibody, MAC-1 Antibody, MO1A Antibody, Neutrophil adherence receptor Antibody, MAC1A Antibody, SLEB6 Antibody, CR-3 alpha chain Antibody
Family / Subfamily:	Integrin / not assigned-Integrin
Host	ITGAM antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	F6.2
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ITGAM / CD11b antibody was raised against Human
Specificity:	Human ITGAM / CD11b
Reactivity:	Human
Purification:	Protein A/G purified
Presentation:	PBS, 0.08% sodium azide, 0.2% carrier protein, sterile-filtered
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	PBMC: Add 10 ul of antibody/10 ⁶ PBMC in 100 ul PBS. Mix gently and incubate for 15 minutes at 2 to 8°C. Wash twice with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. WHOLE BLOOD: Add 10 ul of antibody/100 ul of whole blood. Mix gently and incubate for 15 minutes at room temperature (20°C). Lyse the whole blood. Wash once with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. See instrument manufacturers instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	