

CD45 + CD34 Mouse anti-Human Monoclonal (FITC, PE) (ESS) Antibody - LS-C140428 - LSBio	
CatalogID:	LS-C140428
Target:	CD45/CD34
Host	CD45/CD34 antibody was produced in Mouse
Clonality:	Monoclonal
Clone Name:	ESS
Conjugations:	Fluorescein (FITC), Phycoerythrin (PE)
Immunogen Species:	CD45 + CD34 antibody was raised against Human
Specificity:	Human CD45/CD34
Reactivity:	Human
Purification:	Protein A/G purified
Presentation:	PBS, 0.08% sodium azide, 0.2% carrier protein, sterile-filtered
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	PBMC: Add 10 ul of antibody/10 ⁶ PBMC in 100 ul PBS. Mix gently and incubate for 15 minutes at 2 to 8°C. Wash twice with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze. WHOLE BLOOD: Add 10 ul of antibody/100 ul of whole blood. Mix gently and incubate for 15 minutes at room temperature (20°C). Lyse the whole blood. Wash once with PBS and analyze or fix with 0.5% v/v of paraformaldehyde in PBS and analyze.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst or 100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	