

CD2 + CD19 Mouse anti-Human Monoclonal (FITC, PE) (3.22) Antibody - LS-C140448 - LSBio

CatalogID:	LS-C140448
Target:	CD2/CD19
Host	CD2/CD19 antibody was produced in Mouse
Clonality:	Monoclonal
Clone Name:	3.22
Conjugations:	Fluorescein (FITC), Phycoerythrin (PE)
Immunogen Species:	CD2 + CD19 antibody was raised against Human
Immunogen:	CD2 + CD19 antibody was raised against t Lymphocytes activated by mixed lymphocyte culture.
Specificity:	Human CD19
Reactivity:	Human
Purification:	Protein A/G purified
Presentation:	PBS, 0.08% sodium azide, 0.2% carrier protein, sterile-filtered
Recommended Storage:	Store at 4°C. Do not freeze.
Usage Summary:	PBMC: Add 10 ul of antibody/10 ⁶ PBMC in 100 ul PBS. Mix gently and incubate for 15 minutes at 2 to 8°C. Wash twice with PBS and analyze. WHOLE BLOOD: Add 10 ul of antibody/100 ul of whole blood. Mix gently and incubate for 15 minutes at room temperature 20°C. Lyse the whole blood. Wash once with PBS and analyze. See instrument manufacturers instructions for Lysed Whole Blood and Immunofluorescence analysis with a flow cytometer or microscope.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst or 100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences