

Phosphotyrosine Mouse Monoclonal (FITC) (PY20) Antibody - LS-C140470 - LSBio

CatalogID:	LS-C140470
Target:	Phosphotyrosine
Host	Phosphotyrosine antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	PY20
Conjugations:	Fluorescein (FITC)
Specificity:	Phosphotyrosine
Presentation:	20 mM sodium phosphate, 150 mM sodium chloride, 50% glycerol, 3 mM sodium azide, sterile-filtered
Recommended Storage:	Store at -20°C. Aliquot to avoid freeze/thaw cycles.
Usage Summary:	Application: A dilution of 5-10 ug/ml is suggested for immunohistochemistry or immunocytochemistry.
Uses:	IHC, Western blot (Optimal dilution to be determined by the researcher)
Size:	50 µg or 100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences