

ADM2 / Adrenomedullin 2 Rabbit anti-Rat Polyclonal (Rhodamine) Antibody - LS-C143624 - LSBio	
<b>CatalogID:</b>	LS-C143624
<b>Target:</b>	Adrenomedullin 2 (ADM2)
<b>Synonyms:</b>	ADM2 Antibody, AM2 Antibody, DJ579N16.4 Antibody, Adrenomedullin 2 Antibody, Intermedin Antibody
<b>Host</b>	ADM2 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Rhodamine
<b>Immunogen Species:</b>	ADM2 / Adrenomedullin 2 antibody was raised against Rat
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ADM2 / Adrenomedullin 2 antibody was raised against synthetic peptide corresponding to rat Intermedin
<b>Specificity:</b>	Rat Intermedin.
<b>Reactivity:</b>	Rat
<b>Purification:</b>	Purified
<b>Presentation:</b>	Liquid
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	