

ADM2 / Adrenomedullin 2 Rabbit anti-Human Polyclonal (aa57-92) (Rhodamine) Antibody - LS-C143631 - LSBio

CatalogID:	LS-C143631
Target:	Adrenomedullin 2 (ADM2)
Synonyms:	ADM2 Antibody, AM2 Antibody, DJ579N16.4 Antibody, Adrenomedullin 2 Antibody, Intermedin Antibody
Host	ADM2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Rhodamine
Immunogen Species:	ADM2 / Adrenomedullin 2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ADM2 / Adrenomedullin 2 antibody was raised against synthetic peptide corresponding to aa57-92 of human Intermedin
Specificity:	Human Intermedin.
Epitope:	aa57-92
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences