

Catestatin Rabbit anti-Mouse Polyclonal (FITC) Antibody - LS-C143692 - LSBio

CatalogID:	LS-C143692
Target:	Catestatin
Host	Catestatin antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Catestatin antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	Catestatin antibody was raised against synthetic peptide corresponding to mouse Catestatin
Specificity:	Mouse Catestatin.
Reactivity:	Mouse, Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences