

NERP-1 / Neuroendocrine Regulatory Peptide-1 Rabbit anti-Human Polyclonal (aa281-306) (FITC) Antibody - LS-C144037 - LSBio	
CatalogID:	LS-C144037
Target:	NERP-1 / Neuroendocrine Regulatory Peptide-1
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Human
Antigen Type:	Synthetic peptide
Immunogen:	Synthetic peptide corresponding to aa281-306 of human NERP-1
Specificity:	Human NERP-1.
Epitope:	aa281-306
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014 © 2014 LifeSpan BioSciences	