

NPB29 / Neuropeptide B-29 Rabbit anti-Rat Polyclonal (Rhodamine) Antibody - LS-C144059 - LSBio

CatalogID:	LS-C144059
Target:	NPB29 / Neuropeptide B-29
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Rhodamine
Immunogen Species:	Rat
Antigen Type:	Synthetic peptide
Immunogen:	Synthetic peptide corresponding to rat Neuropeptide B-29
Specificity:	Rat Neuropeptide B-29.
Reactivity:	Rat
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences