

Salusin-Alpha Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C144164 - LSBio

CatalogID:	LS-C144164
Target:	Salusin-Alpha
Host	Salusin-Alpha antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Salusin-Alpha antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Salusin-Alpha antibody was raised against synthetic peptide corresponding to human Salusin-alpha
Specificity:	Human Salusin-alpha.
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences