

SYT4 Rabbit anti-Human Polyclonal (aa73-97) (Cy5) Antibody - LS-C144212 - LSBio

CatalogID:	LS-C144212
Target:	synaptotagmin IV (SYT4)
Synonyms:	SYT4 Antibody, KIAA1342 Antibody, Synaptotagmin IV Antibody, Synaptotagmin-4 Antibody, SytIV Antibody, HsT1192 Antibody, Synaptotagmin 4 Antibody
Host	SYT4 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Cy5
Immunogen Species:	SYT4 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SYT4 antibody was raised against synthetic peptide corresponding to aa73-97 of human Synaptotagmin IV
Specificity:	Human Synaptotagmin IV.
Epitope:	aa73-97
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences