

ULBP3 Mouse anti-Human Monoclonal (PE) Antibody - LS-C144249 - LSBio	
CatalogID:	LS-C144249
Target:	UL16 binding protein 3 (ULBP3)
Synonyms:	ULBP3 Antibody, ALCAN-gamma Antibody, NKG2D ligand 3 Antibody, NKG2DL3 Antibody, RAET1N Antibody, N2DL3 Antibody, UL16 binding protein 3 Antibody, UL16-binding protein 3 Antibody, N2DL-3 Antibody
Host	ULBP3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ULBP3 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	ULBP3 antibody was raised against baF3 cells transfected with human ULBP-3 (hULBP-3)
Specificity:	Detects human ULBP-3.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	Saline, up to 0.5% BSA, 0.1% sodium azide
Recommended Storage:	May be stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Dilute only prior to immediate use. Stable for at least 12 months at 4°C. Sensitive to light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	