

UCN3 / SPC Rabbit anti-Human Polyclonal (Cy5) Antibody - LS-C144260 - LSBio	
<b>CatalogID:</b>	LS-C144260
<b>Target:</b>	urocortin 3 (UCN3)
<b>Synonyms:</b>	UCN3 Antibody, SPC Antibody, Stresscopin Antibody, Urocortin 3 (stresscopin) Antibody, Urocortin 3 Antibody, Urocortin III Antibody, Urocortin-3 Antibody, UCNIII Antibody, SCP Antibody, Ucn III Antibody
<b>Host</b>	UCN3 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Cy5
<b>Immunogen Species:</b>	UCN3 / SPC antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	UCN3 / SPC antibody was raised against synthetic peptide corresponding to human Urocortin III
<b>Specificity:</b>	Human Urocortin III.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	Liquid
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/25/2014	
© 2014 LifeSpan BioSciences	