

UCN3 / SPC Rabbit anti-Mouse Polyclonal (Rhodamine) Antibody - LS-C144263 - LSBio

CatalogID:	LS-C144263
Target:	urocortin 3 (UCN3)
Synonyms:	UCN3 Antibody, SPC Antibody, Stresscopin Antibody, Urocortin 3 (stresscopin) Antibody, Urocortin 3 Antibody, Urocortin III Antibody, Urocortin-3 Antibody, UCNIII Antibody, SCP Antibody, Ucn III Antibody
Host	UCN3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Rhodamine
Immunogen Species:	UCN3 / SPC antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	UCN3 / SPC antibody was raised against synthetic peptide corresponding to mouse Urocortin III
Specificity:	Mouse Urocortin III.
Reactivity:	Mouse
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences